# Clenny Aerosol<sup>®</sup>


PRODOTTO DA:  
MANUFACTURED BY:

CE 0123

**SISTEMA PER AEROSOLTERAPIA A MICRO COMPRESSORE  
MICRO COMPRESSOR NEBULIZER SYSTEM**

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## EFFICIENCY OF THE BREATH-ENHANCED NEBULIZER

The expiration time is considered to be 60% of the complete respiratory act. As a result of this, in all "Continuous Flow" systems **1** more than 60% of drug is wasted in the air.

The exclusive CLENNY AEROSOL double-valve system reduces the wasted medication, increases the **2** output and adjusts the delivery of drug according to the patient breathing.

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## FAST DRUG DELIVERY

The wide inspiration channel on one side **1** draws extra air from outside by Venturi effect **2**, and together with air flow generated by the compressor **3**, the quality and fastness of the drug output increase significantly **4**.

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
## MORE RESPIRABLE FLOW

**1 Inspiration Phase:** the inspiration valve opens drawing outside air by Venturi effect and increases the nebulizing capacity. The expiration valve is closed to create the necessary depression to open the inspiration valve.

**2 Expiration Phase:** the inspiration valve closes reducing the nebulizing capacity. The expiration valve opens for an easy expiration.

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
## QUALITY IN NEBULIZATION

The efficiency of the nebulizer and the power of the micro compressor deliver particles below 5µm to reach the very deep airways and increase the respirable volume. CLENNY AEROSOL nebulizes a saline solution, NaCl 0,9%, with MMAD 1,9 µm (tested with Api Aerosizer Mach 2) or MMD 4,7 µm (Malvern Spraytec/Insittec). CLENNY AEROSOL was tested with the most popular drugs.

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## 75% OF ACTIVE PRINCIPLE RESPIRED DURING THE FIRST 3 MINUTES

A test conducted in Rome, by "La Sapienza" University, compared CLENNY AEROSOL with a traditional Continuous Flow system, measuring with a filter the quantity of solution respired. This test showed a higher deposition with CLENNY AEROSOL, and 75% of the active principle deposited at the 3rd minute, while the plateau is achieved at the 5th minute.

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Data achieved with Beclometasone 0,04%.

**MEASURING METHOD:** data have been obtained with a breathing simulator (16 l/min for 15 breaths) and a filter in the inhalation phase able to retain 99,99% of the nebulized material. The quantity of active principle deposited onto the same filter has been determined with HPLC test.

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## STANDARD ACCESSORIES AND OPTIONALS

The patient compliance towards aerosol therapy depends on the convenience and simplicity of the device. At the same time, the accessories supplied with the unit must comply with the therapeutical requirements. All CLENNY AEROSOL accessories have been designed to increase the patient compliance respecting the safety and hygiene characteristics necessary for a medical device.

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## MICROCOMPRESSOR NEBULIZER SYSTEM

- *Small and convenient device, to be used in any situation.*
- *Powerful and silent microcompressor.*
- *Rapid and efficient nebulizer, thanks to the new inhalation/exhalation valve system which increases the medication administration.*
- *Mouthpiece storable in the nebulizer (2 pcs. included).*
- *Non-invasive nosepiece.*
- *State-of-the-art adult and pediatric masks for the proper operation of the valve system even with non-cooperating patients.*
- *Air filter.*
- *Air tubing link between the nebulizer and the microcompressor.*

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MEDEL S.p.A. • VIA MICHELI 7/9, 43056 S. POLO DI TORRILE (PARMA) ITALY