



AC15

ASPIRATOR FOR POST TRACHEOTOMY

- Ideal for home post-surgery treatments.
- Reliable, maintenance free and silent MEDEL pump.
- Shock-resistant, easy-to-carry and easy-to-clean chassis.
- Suction adjustment and vacuum gauge.
- 1,200 ml. polycarbonate reusable canister with overflow valve.
- Antibacterial filter.



Technical Specificaton

Size:	350 x 180 x 200 mm
Weight:	3 Kg
Current:	1 A
Max vacuum:	-82 kPa; -610 mmHg
Max flow:	15 l/min
Voltage:	230V, 50Hz. Other voltages available upon request
Warranty:	2 years
Manufacturer:	Medel S.p.A

Ref: 90746 - CE0123

Technical specifications and appearance of products are subject to change without notice.