



DC15

ASPIRATOR FOR AMBULANCE / BATTERY OPERATED

Ideal for ambulance and homecare.

Three versions:

Model DC15 BASIC: 12VDC, supplied with car or ambulance cable only.

Model DC15 UNIVERSAL: 12VDC, supplied with car or ambulance cable and universal voltage (100 – 240VAC) adapter.

Model DC15 BATTERY: 12VDC, supplied with car or ambulance cable, universal voltage (100 – 240VAC) adapter and internal rechargeable battery pack, 45 minutes run time.

- Reliable, maintenance free and silent MEDEL pump.
- Shock-resistant, easy-to-carry and easy-to-clean chassis.
- Suction adjustment and vacuum gauge.
- 1,200 ml. polycarbonate reusable canister with overflow valve.
- Antibacterial filter.



Technical Specificaton

Size:	350 x 180 x 200 mm
Weight:	2,6 Kg(3,2 kg with battery)
Current:	2,5 A
Max vacuum:	-82 kPa; -610 mmHg
Max flow:	15 l/min
Battery*:	3300mAh(Ni-MH battery cells) *only for DC15 battery
Warranty:	2 years
Manufacturer:	Medel S.p.A
Ref: 90823 - DC15 Basic - CE0123	
Ref: 90832 - DC15 - CE0123	
Ref: 90824 - DC15 Battery - CE0123	