

medel[®]

DAILY

...for your health



PULMODATA

d.o.o.

Profesionalna medicinska oprema
razvoj, inženiring, zastopanje

Balantičeva 9
1234 MENGEŠ, SLOVENIJA
tel.: 01/723 95 11, fax: 01/723 95 15
e-mail: info@pulmodata.si
www.pulmodata.si

CE 0123

**MISURATORE DI PRESSIONE AUTOMATICO DA POLSO
AUTOMATIC WRIST WATCH BLOOD PRESSURE MONITOR**



DAILY

...for your health



DISPLAY BOX:
6 PEZZI/PCS.

MISURATORE DI PRESSIONE AUTOMATICO DA POLSO

La tecnologia Fuzzy Logic consente di eliminare le operazioni di impostazione del misuratore. Infatti, già in fase di gonfiaggio viene effettuata la misurazione della pressione sistolica e automaticamente viene regolato il livello della pressione di gonfiaggio del bracciale. Innovativo sistema di sgonfiaggio "lineare" con decremento costante e controllo elettronico: garantisce una misurazione più attendibile e precisa

- Misurazione della pressione sistolica, diastolica e delle pulsazioni
- Precisione clinicamente validata
- Memoria delle ultime 30 misurazioni
- Registrazione di data ed ora associate a ciascuna misurazione
- Gonfiaggio e sgonfiaggio automatico del bracciale
- Coperchio antiurto a salvaguardia del display e di accensioni accidentali ed involontarie
- Grande display a cristalli liquidi per una agevole lettura
- Comoda borsina di trasporto

**GARANZIA 5 ANNI
5 YEARS WARRANTY**

AUTOMATIC WRIST WATCH BLOOD PRESSURE MONITOR

The Fuzzy Logic technology allows eliminating pre-selection operations. In fact, as early as during inflating phase, the blood pressure monitor measures the systolic pressure and regulates automatically the cuff inflating pressure level.

New "linear" deflation system with constant deflating phase and electronic check for an accurate and reliable measurement.

- Systolic, diastolic and pulse rate measurements
- Clinically validated
- 30 readings memory capacity
- Fully automatic inflation and deflation
- Active flip cover turns off the device automatically to extend battery life and to protect LCD
- Large easy-to-read display
- Stylish and convenient carrying case

Caratteristiche Tecniche - Technical Specifications

MEDEL DAILY

Dimensioni - Size:	84 mm x 75 mm x 70 mm
Peso - Weight:	150 g
Intervallo pressione - Range SYS/DIA:	30-280 mmHg
Pulsazioni - Pulse:	40-200/min.
Precisione pressione - Pressure accuracy:	+/- 3 mmHg
Alimentazione - Power source:	2 batterie UM-3, dim. AA, 1.5V • 2 batteries UM-3, dim. AA, 1.5V
Garanzia - Warranty:	5 anni • 5 years

REF 90485 - MEDEL DAILY

DISPLAY BOX: n° 6 pz • no. 6 pcs.

Dimensioni • Size: 220x270x100 mm

Peso • Weight: 2 kg

MASTER: n° 6 display box • no. 6 display boxes

n° 36 pz • no. 36 pcs.

Dimensioni • Size: 470x300x350 mm



8 025081 904858

Conforme ai requisiti della direttiva 93/42/CEE Dispositivo Medico Classe IIa • Complies with 93/42/EEC Directive. Class IIa Medical Device

MEDEL S.p.A. • VIA MICHELI 7/9, 43056 S. POLO DI TORRILE (PARMA) ITALY
TEL.(+39) 0521 315211 - FAX (+39) 0521 315212 - E-MAIL: INFO@MEDEL.IT - WWW.MEDEL.IT