



FAMILY

...for your health



CE 0123

**SISTEMA PER AEROSOLTERAPIA PER USO DOMICILIARE
COMPRESSOR NEBULIZER SYSTEM FOR HOME USE**



FAMILY

...for your health

SILVER



2 BLUE



MEDEL FAMILY è un sistema ultra-compatto per aerosolterapia ad aria compressa. La struttura di ridotte dimensioni e la praticità d'uso fanno di MEDEL FAMILY un apparecchio ideale per l'utilizzo domiciliare.

MEDEL FAMILY è dotato del nebulizzatore MEDELJET BASIC che garantisce una maggior rapidità ed efficacia di nebulizzazione. La versione MEDEL FAMILY 2BLUE è fornita in 2 gradevoli colori mentre la versione MEDEL FAMILY SILVER è completa di una comoda borsa di trasporto.

MEDEL FAMILY è dotato di: nebulizzatore MEDELJET BASIC, boccaglio, raccordo a T, forcilla nasale, mascherina adulti, tubo aria, cavo alimentazione, 5 filtri aria.

MEDEL FAMILY is an ultra-compact compressor nebulizer system. The compact structure and ease to use make MEDEL FAMILY ideal for home use.

MEDEL FAMILY is equipped with MEDELJET BASIC nebulizer, that guarantees high nebulization rate and efficiency. MEDEL FAMILY 2BLUE version is proposed in two nice colours, MEDEL FAMILY SILVER is equipped with a convenient carrying case.

MEDEL FAMILY is equipped with a complete set of accessories: MEDELJET BASIC nebulizer, mouthpiece, T piece, nosepiece, adult mask, air tube, power cord, 5 spare filters.

Caratteristiche Tecniche - Technical Specifications

MEDEL FAMILY SILVER - 2BLUE

Nebulizzazione: - Nebulization rate:

ca. 0,28 ml/ min.

Particelle - Particle size:

MMAD 1,9 µm (NaCl 0,9%) - API AEROSIZER
MMD 4,99 µm - MALVERN SPRAYTEC/INSITEC

Pressione max. al compressore - Compressor max. pressure:

2,30 bar

Peso - Weight:

1,55 kg

Dimensioni - Size:

167 mm x 145 mm x 97 mm

Livello sonoro - Noise level:

55 dBA

REF 90543 - MEDEL FAMILY SILVER

REF 90479 - MEDEL FAMILY 2BLUE

MASTER: n° 6 pz • no. 6 pcs.

Dimensioni • Size: 740 x 230 x 230(h) mm

Peso • Weight: 12,7 kg

PALLET: n° 30 master • no. 30 masters

n° 180 pz • no. 180 pcs.



MASTER: n° 6 pz • no. 6 pcs.

Dimensioni • Size: 470 x 220 x 290(h) mm

Peso • Weight: 11,7 kg

PALLET: n° 30 master • no. 30 masters

n° 180 pz • no. 180 pcs.



Conforme ai requisiti della direttiva 93/42/CEE Dispositivo Medico Classe IIa • Complies with 93/42/EEC Directive. Class IIa Medical Device

MEDEL S.p.A. • VIA MICHELI 7/9, 43056 S. POLO DI TORRILE (PARMA) ITALY
TEL.(+39) 0521 315211 - FAX (+39) 0521 315212 - E-MAIL: INFO@MEDEL.IT - WWW.MEDEL.IT