

# Biocomfort



**massator premium**

**Massage using the  
massator premium for  
overall well-being**

# biocomfort massator premium



- Attractive design
- Ergonomically shaped handle
- Convenient coiled cable
- Easy to grip surface



## Unique:

The infinitely adjustable electronic control provides a smooth transition from gentle to strong massage

Simple and easy to use switch and controller for convenient handling



## Scope of delivery/available accessories

- **Base part of storage case** with insert for foot massage
- **Cover of storage case** for accessories
- **Massator premium** with heatable base plate for body massage
- **Natural bristle plate** for body massage
- **Massage brush** for scalp, arms and legs
- **Nub plate** for thighs and soles of the feet
- **Cosmetic brush** for face and neck
- **Cosmetic sponge** for face and neck
- **Suction cup** for face, neck and bust
- **Massage finger** for reflex-zone massage and acupuncture
- **Base plate with Velcro® covering** for Shiatsu attachments
- **Shiatsu set** for Shiatsu massage
- **Joint-massage suction cup** particularly for massage of joints and body massage
- **Stand** for independent back massage without the need for additional help
- **Wall bracket** for stand including screws and dowels

## Technical data:

- Mains voltage: 230 V, 50 Hz
- Motor power: 45 W
- Motor speed: from approx. 700 to 2750 rpm
- Heating plate power: 50 W
- Preheating of massage plate from approx. 20° to 45°C approx. 50 sec.
- Temperature of base plate with activated heating: Upper temperature range approx. 45°C Lower temperature range approx. 35°C
- Effective cable length: 3 meter
- Weight of massator premium: 1,4 kg
- Weight of storage case: 1,5 kg
- Height adjustment range on stand: 290 mm
- Angle adjustment range on stand: 60°



For convenient use while sitting in an easy chair or any other seating accommodation. The massage will help to treat cramps and hardening of the muscles. For overall well-being.



For the back massage we recommend the use of the massage brush, the natural bristle brush or the joint-massage suction cup – whichever suits you best.



Effective and convenient reflex-zone massage for the soles of your feet using the screwed-on massage finger.

## Your contact person